DISCOCERV™


Cervical Disc Prosthesis

SURGICAL TECHNIQUE GUIDE


Scient'x[®] *α* Alphatec Spine[®]
Solutions for the Aging Spine[®]

Discocerv™: the evolution of the Cervidisc cervical disc prosthesis

Operating time saving

Color coding to optimize instruments identification according to the size of the prosthesis

Optimized design

Convex upper and lower plates to comply with the anatomical curvature of the disc space

Longevity

The plates are made from **titanium** TA6V ELI (ASTM F136 or ISO 5832-3) which gives the prosthesis excellent mechanical and biomechanical properties

The superior ceramic part (convex) is made from **alumina** (AL2O3 ISO 6474 type A) and the inferior part (concave) from **zirconia** (ZrO₂HfO₂Y₂O₃ ISO 13356) giving excellent wear resistance in the long term

The cervical disc prosthesis DISCOCERV™ – Cervidisc Evolution is a surgical device intended to preserve the physiological mobility of the cervical functional spinal unit via the anterior route.

The disc replacement is only indicated for levels from C3 to C7. The main indication is disc degeneration without instability.

The indications include the following :

- degenerative disc disease
- herniated disc
- myelopathy associated with a spondylotic stenosis of the foramen or canal
- root disease associated with a neurological deficiency that does not respond to conservative treatments

Mobility

Angulation of 9° in the sagittal and coronal planes to preserve the physiological amplitude of a normal disc


Security

Bi-directional grooves to avoid antero-posterior and lateral migration of the titanium plates

Controlled movement

Perfect joining of the inferior and superior ceramic parts so as to facilitate prosthesis mobility

Pre-operative considerations

Indications :

The DISCOCERV™ cervical disc prosthesis is a surgical implant intended to maintain the physiological cervical mobility of a spinal segment. The prosthesis is implanted between C3 and C7 by the anterior route.

The principal indication is degeneration of the disc without instability. The indications are as follow :

- degenerative disc disease
- herniated disc
- myelopathy associated with a spondylotic stenosis of the foramen or canal
- root disease associated with a neurological deficiency that does not respond to medical treatments

Contraindications :

The contraindications to the cervical disc prosthesis DISCOCERV™ are as follow :

- active infection or local inflammation
- vertebral osteoporosis
- vertebral tumor disease
- injury
- local deformity
- instability
- surgical history at the site to be treated
- rheumatoid arthritis
- metabolic bone disease
- pregnancy
- incompatible age or physical condition of the patient
- obvious lack of mobility at the level concerned
- allergy or intolerance towards the materials constituting the device
- any situation not included in the indications

The DISCOCERV™ cervical disc prosthesis must not be used in the context of disorders requiring a multi-staged arthroplasty or an adjacent arthrodesis (above or below).

The DISCOCERV™ cervical disc prosthesis is not designed, intended or sold for uses other than those indicated.



Pre-operative X-rays

A pre-operative scanner section and profile X-ray should preferably be performed to determine the ideal size of DISCOCERV™ cervical disc prosthesis to be implanted.



1 Surgical protocol/patient positioning

- a** The surgery should be performed under general anesthesia. The patient should be placed in dorsal decubitus on an ordinary table. The head should be placed in the neutral position (no flexion or hyperextension) ; it may be held in place with an adhesive band. It is vital to keep a strict profile to avoid losing the median line. **The success of surgery depends on this.** It is possible to position the head on a transparent headrest or a transparent extension of the table.
- b** The approach is a standard right anterior pre-sterno-cleido-mastoid approach (Cloward or Smith-Robinson).
- c** Under scopic control, determine the correct intervertebral level (image intensifier).
- d** After incision of the anterior longitudinal ligament, perform the disc resection using appropriate rongeurs and/or curettes.

2 Discectomy and preparation of the surgical site

When the disc has been excised, proceed with a bilateral postero-lateral foraminal decompression. If the disc space is too tight, an intersomatic retractor may be used. Decompression may require a high speed drill or a kerisson to free the bone (osteophytosis). This surgical gesture should be performed under surgical microscope.

3 Preparation of the endplates

The vertebral plates should be prepared and cleaned carefully with the curette. A damaged cartilage may cause subsidence of the prosthesis in the vertebral body. This procedure is essential to ensure optimum contact between the implant and adjacent vertebral plates, to avoid secondary displacement.

Note : It is important to insert the two flat probes one after the other to determine the size of the most suitable prosthesis. In order to support a better distribution of the constraints on the vertebral endplates, the prosthesis covering the largest part of the endplate should be chosen.

4 Control of the intervertebral space

Once the site preparation is finalized, the **flat probe** should be inserted in the intervertebral space. The flat probe allows to verify the quality of the intervertebral space created after discectomy and select the correct prosthesis footprint.

The flat probes feature a midline laser marking on its superior part that can help to evaluate the midline of the intervertebral space.

Based on this midline laser marking; a marking realized with help of drilling tool, onto the superior and inferior vertebral body can help to keep the midline till the end of the surgery.

The flat probes also feature a color coding helping to identify the prosthesis footprint (Yellow: Small footprint; Purple: Large footprint)

Note : The trials can be introduced and removed easily since their surface have no retaining edges. Trials reflect the height of the implants without teeth. Trials as flat probes also feature laser marking to identify midline of the intervertebral space

5 Insertion of the prosthesis trials

After selection of the appropriate footprint, the set of **prosthesis trials** must be used to determine the ideal size of the prosthesis to be inserted. The trials also feature the same color coding system as the flat probes, according to the footprint selection in previous step, only the trials following the same color coding will be used. The design of the trials perfectly matches the dimension of the prosthesis without anchoring teeth. The trials are used to determine the implant best suited to the patient's anatomy and requirements.

The choice of the right size of the prosthesis is the most critical part. If the prosthesis is too small, it may migrate and if it is too big, it may interfere with mobility.

Flat probe
21GPC0_

Prosthesis trials
21FPDC_-H_

To determine the size of the most suitable implant, it is recommended to first insert the smallest trial implant and sequentially select the next biggest one till the proper size is achieved.

The height to be restored can be determined according to the adjacent disc heights or the facet parallelism at the index level. This is left to the surgeon's discretion.

Following the reference and laser marking indicated onto the trial, the corresponding appropriate prosthesis is selected.


6 Insertion of the cervical distractor

The **distractor** should then be inserted in the intervertebral space. It is recommended to slightly distract the intervertebral space during the prosthesis implantation time. If the anatomy of the patient does not allow the introduction of the specific distractor, a Caspar distractor can be used.

The distractor should be correctly positioned in the transversal plane and should be centered along the midline. The distractor branches should be in contact with the posterior wall.

Scopic control should be performed to ensure correct placement of the distractor:

- frontal control is used to verify correct positioning of the distractor in the transversal plane and exact centering along the midline.
- profile control is used to verify the retraction of the vertebral plates and if needed the retraction of the joint facets.

Note : To optimize the surgical gesture the middle line may be marked on the lower and upper vertebral plates.


Distractor
21EDCO_


Holder
21PRE0_

TOP VIEW OF THE HOLDER

FRONTAL VIEW OF THE HOLDER


Note : the inscription "TOP" should always be on the top of the instrument

7 Positioning of the prosthesis on the holder

The **holder** is used for impaction of the prosthesis in the intervertebral space. The holder tips fit on the lateral sides of the prosthesis to ensure a firm grasp of the implant.

Firstly, the upper plate (including spherical ball) of the prosthesis is positioned manually (A) on the two upper arms of the holder. The lower plate of the prosthesis is then inserted (B) in the two lower arms of the holder using the flexibility of the instrument.

Note : it is important to obtain a perfect coupling between the two ceramic parts.

Note : To ensure a correct positioning of the implant onto the holder, make sure the frontal part of the prosthesis (including laser marking logo and references) are in contact with the holder. When inserted into the patient, these laser marking logo and references should be visible by the surgeon.

Note : No lever or torsion movement should be applied to the holder during impaction of the prosthesis.

8 Insertion of the prosthesis

The holder takes up little space because it is in line with the prosthesis, leaving the operative zone totally visible. The distractor is used to guide the holder and facilitate visualization of the correct alignment of the prosthesis. The effort during insertion of the prosthesis should be along the ancillary axis.

Superior key
21CLE03 and
21CLE05

Inferior key
21CLE04


The posterior edge of the prosthesis must be positioned as close as possible to the posterior wall of the operated vertebrae. That is why it is necessary to correctly prepare the intervertebral space and to remove all osteophytes.

When the prosthesis is in place, the distractor and then the holder may be withdrawn. It is important to check that the prosthesis is correctly placed with a per-operative X-ray.

A proper implant positioning means that the prosthesis is perfectly aligned with the posterior wall.

The inferior and superior key can help to adjust the final positioning of the implant.

Note : After scopic control and if necessary, the superior and inferior keys may be used to reposition the prosthesis. For clear identification of the repositioning tool concerning the footprint size, the superior key features the same color coding as previously described (yellow: small footprint; purple: bigger footprint).

9 Closure of the approach wound

The approach wound is closed after rinsing. Hemostasis is verified and drainage may be required. The platysma should be carefully reconstituted and the skin closed by absorbable intradermic sutures.

Removing the implant

If it becomes necessary to remove the prosthesis, the cervical approach should be performed again until exposure of the instrumented zone. The distractor should be inserted and the **two retrieval keys** placed on the prosthesis to facilitate withdrawal of the implant.

Retrieval key
21PIN04

DISCOCERVTM

Cervidisc evolution

INSTRUCTIONS, IMPLANTS & INSTRUMENTS

Instruction for use

■ OBJECTIVE

The Cervical Disc Prosthesis DISCOCERV is an implant intended for the surgical treatment of the cervical spine. The Cervical Disc Prosthesis DISCOCERV is intended for patients suffering from degenerative disorders of the cervical spine, but who do not display any instability.

■ GENERAL DESCRIPTION

The Cervical Disc Prosthesis DISCOCERV has been designed in order to preserve cervical mobility. The preserved mobility is similar to physiological mobility. The Cervical Disc Prosthesis DISCOCERV takes the form of a ceramic joint, consisting of a cup and a spherical ball. The two parts that form the joint are embedded in metal plates that come into contact with the vertebral bodies. The upper surface of the disc replacement inserted into the disc chamber is convex in the sagittal plane, and has a convex profile on its lower surface, in order to match the vertebral plate to the frontal plane as closely as possible. Several different thicknesses are supplied in order to fit different morphologies of the intervertebral spaces. The Cervical Disc Prosthesis DISCOCERV is manufactured from ELI titanium (according to ISO5832-3 or ASTM F136), and from two different types of ceramics: zirconia and alumina (ISO 13356 and ISO 6474, respectively). The titanium endplates have been plasma sprayed (with Titanium according to 5832-3 or ASTM F67).

The Cervical Disc Prosthesis DISCOCERV must not be used with components supplied by other manufacturers.

■ INDICATIONS

The Cervical Disc Prosthesis DISCOCERV is a surgical implant intended to preserve the physiological cervical mobility of a segment of the spine via the anterior route. The disc replacement is only indicated for levels from C3 to C7. The main indication is disc degeneration that does not display any instability.

The indications include the following:

- degenerative disc disease,
- herniated disc,
- myelopathy associated with a spondylotic stenosis of the foramen or canal,
- root disease associated with a neurological deficiency that does not respond to medical treatments.

■ CONTRAINDICATIONS

The contraindications for the Cervical Disc Prosthesis DISCOCERV are as follows:

- active infection or local inflammation,
- vertebral osteoporosis,
- vertebral tumor disease,
- injury,
- local deformity,
- instability,
- surgical history at the site to be treated,
- rheumatoid arthritis,
- metabolic bone disease,
- pregnancy,
- incompatible age or physical condition of the patient,
- obvious lack of mobility at the level concerned,
- allergy or intolerance towards the materials constituting the device,
- any situation not included in the indications.

The Cervical Disc Prosthesis DISCOCERV must not be used in the context of disorders requiring a multi-staged arthroplasty or an adjacent arthrodesis (above or below).

The Cervical Disc Prosthesis DISCOCERV is not designed, intended or sold for any uses other than those indicated.

■ POSSIBLE SIDE EFFECTS

- infection,
- post-operative migration of the implant,
- spontaneous fusion,
- in cases of fusion affecting the vertebral stages adjacent to the level of the arthroplasty,
- dislocation,
- the plates of the disc replacement may be driven down into the vertebral bodies,
- intolerance towards the material.

Note: An additional surgical operation may be necessary to correct a side effect.

Warnings: Not every surgical operation leads to an entirely satisfactory outcome. This is particularly true of spinal surgery, where numerous external factors can compromise the results.

■ SURGICAL PRECAUTIONS

The surgeon must be fully conversant with the Cervical Disc Prosthesis DISCOCERV, the method of application, the instruments and the surgical procedure. The Cervical Disc Prosthesis DISCOCERV must be implanted using the recommended operative technique. The appropriate size of the Cervical Disc Prosthesis DISCOCERV must be selected according to the height of the disc to be restored, a height that is inappropriate for the height of the normal interbody space could compromise the clinical outcome. It is recommended that an intervertebral distractor be used to carry out the discectomy and insertion of the disc replacement. Before

implanting the Cervical Disc Prosthesis DISCOCERV, the vertebral plates must be carefully prepared by removing the cartilaginous layer and roughening the subchondral bone. It is important not to remove any of the bony endplate as this could lead to subsidence of the Cervical Disc Prosthesis DISCOCERV. Obtaining an absolutely stable primary anchorage of the disc replacement is essential in order to achieve lasting fixation of the implant. The following defects can lead to dislodgement, subsidence and compression of the implant or other post-operative complications:

- excessive removal of bone when preparing the endplate,
- selection of an inappropriate size of disc prosthesis,
- forcible introduction of the replacement, resulting in vertebral fractures,
- be careful to avoid excessive resection of the uncus, which could weaken the vertebral plate.

Secondary anchorage is achieved by means of osteointegration that takes place within the porous structure of the titanium plates of the replacement. This porous structure is obtained by applying a deposit of titanium T40 (ISO 5832-3 or ASTM F67). An X-ray checkup makes it possible to confirm that the cervical disc replacement has been correctly positioned relative to the vertebrae. After the implantation, the batch number and the reference code of the Cervical Disc Prosthesis DISCOCERV implanted must always be recorded in the patient's surgical file. This can be done using the labels provided inside the sterile packaging.

This product is a single use device. Under no circumstances should it be reused. While the device may appear to be undamaged, it may have small defects or internal stress patterns, as a result of the prior implantation or removal that could lead to fatigue failure. Additionally, please note that the removed implant has not been designed or validated so as to allow for decontamination of microorganisms. Reuse of this product could lead to cross-infection and/or material degradation as a result of the decontamination process. The company accepts no responsibility for products which have been reused.

Attention: no silicone, formaldehyde or latex based products should be implanted with a SCIENT'X implant.

■ PACKAGING

Cervical Disc Prosthesis DISCOCERV are supplied in individual sterile packaging. The associated accessory instruments are supplied in a single container, in kit form. The packaging and the kit must be closed, sealed and intact when received. If the depot kit system is used, it is essential to check that it is complete. The condition of all the implants and instruments must always be checked before use. If the kit is to be reassembled in the depot, the implants are supplied in sterile packs, and the instruments in individual sachets. In this case, the labeling shows the reference and batch number of the product, as well as the details of the manufacturer. Instruments supplied in sachets are not sterile. The storage conditions must be such as to maintain the integrity of the implants, the associated accessories and their respective packaging. If the packaging ensuring the sterility of the implants has been compromised, or if the implants or accessories have been damaged, the devices concerned must not be used, and should be returned to SCIENT'X.

■ DECONTAMINATION, CLEANING AND STERILIZATION

Products delivered in sachets are not sterile.

For implants delivered sterile: the implants are sterilized by Gamma radiation at doses of 25 to 40 kGy. The expiry date is 5 years. The expiry date of sterile parts is indicated on the packaging. Re-sterilization of implants delivered sterile is prohibited.

For instruments not sterile on delivery: all instruments delivered non-sterile must be decontaminated, cleaned and sterilized before and after use. Implants and instruments in sachets must be removed from the original packaging for the following operations:

Recommended method:

- **Decontamination:** Plunge the instruments into a bactericidal and fungicidal solution of the didecylidimethylammonium chloride type diluted to 0.5 % (5mL to 1 litre water). Length of soaking: 20 min. Rinse with demineralized water.
- **Cleaning:** Wash the instruments in a LANCER type machine with suitable cleaning products, rinse and dry. Any product that might damage the equipment is forbidden (such as bleach, formol, etc.).
- **Sterilization:** We recommend sterilising the instruments in an autoclave:
 - pre-heating for 25' at 110°C (1 bar)
 - vacuum 5' (0.8 bar under atmospheric pressure)
 - heating 5' at 120°C (1 bar)
 - vacuum 5' (0.8 bar)
 - sterilization 18' at 134°C (2 bars)
 - drying 20' return to room temperature

■ COMPLAINTS

Any complaints, together with the reference and lot number of the incriminated product, should be sent to Scient'x.

■ ADDITIONAL INFORMATION

For further information, contact SCIENT'X, Only for information Subject to modification without notification For updated IFU, please contact your local customer service.

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
Only for information. Subject to modification without notification. For updated IFU, please contact your local customer service.

Date of revision : 06/2010

> ALL IMPLANTS ARE DELIVERED **STERILE**.

Discocerv™ prosthesis

Height (mm)	Length (mm)	Depth (mm)	Reference
5.25	17	13	11PDC1-H1S
6			11PDC1-H2S
6.75			11PDC1-H3S
7.5			11PDC1-H4S
6	20	15	11PDC2-H1S
6.75			11PDC2-H2S
7.5			11PDC2-H3S
8.25			11PDC2-H4S

Trials for Discocerv™ prosthesis

	Height (mm)	Reference
11PDC1-H1S	5.25	21FPDC1-H1
11PDC1-H2S	6	21FPDC1-H2
11PDC1-H3S	6.75	21FPDC1-H3
11PDC1-H4S	7.5	21FPDC1-H4
11PDC2-H1S	6	21FPDC2-H1
11PDC2-H2S	6.75	21FPDC2-H2
11PDC2-H3S	7.5	21FPDC2-H3
11PDC2-H4S	8.25	21FPDC2-H4

Flat probe for 11PDC1-H_ 21GPC01
Flat probe for 11PDC2-H_ 21GPC02

Distractor for 11PDC1-H_ 21EDC02
Distractor for 11PDC2-H_ 21EDC03

Holder for 11PDC1-H_ 21PRE08
Holder for 11PDC2-H_ 21PRE09

Inferior key 21CLE04

Superior key for 11PDC1-H_ 21CLE03
Superior key for 11PDC2-H_ 21CLE05

Retrieval key 21PIN04

Corporate Headquarters:

ALPHATEC SPINE, INC.
5818 El Camino Real
Carlsbad, CA 92008 - USA
TL +1 460 431 9286
www.alphatecspine.com

International Headquarters:

SCIENT'X-ALPHATEC SPINE
22 avenue Jean Bart
78960 VOISINS-LE-BRETONNEUX - FRANCE
TL +33 (0)1 39 30 69 30 - FX +33 (0)1 30 43 82 77
www.scientx.com - www.alphatecspine.com

This product is not available for sale in the United States of America
Non contractual pictures.
Please read carefully the instructions for use bulletin.
Devices may be subject to modification. Patented.

Scient'x[®] α Alphatec Spine[®]
Solutions for the Aging Spine[®]

Discocerv and the Scient'x-Alphatec Spine logo are registered trademarks of Alphatec Spine. © 2010 Alphatec Spine, Inc. All rights reserved.

DOMKT504 revision 1.1/ DMQ #: 11-421/Approval date: August 2011