





DLIF Cage System




LnK DLIF Cage System

The LnK DLIF Cage can make solid inter-body fusion with optimized anatomical position. Broad bone graft area also makes perfect fusion. It can be used with Neuro-monitoring System to avoid neural damage.



Features & Benefits

- 1) Wide Bone Graft Space
- 2) Variety of angles, lengths, and heights
- 3) Anatomical shape designed
- 4) Radiopaque markers allows for visualization in radiographic image



Implant Specification

Part No.	Dimension (mm)				Part No.	Dimension (mm)				Part No.	Dimension (mm)				Part No.	Dimension (mm)			
	W(M-L)	L(A-P)	A(°)	H(mm)		W(M-L)	L(A-P)	A(°)	H(mm)		W(M-L)	L(A-P)	A(°)	H(mm)		W(M-L)	L(A-P)	A(°)	H(mm)
3810-4008	18	40	0	8	3816-4008	18	40	6	8	3830-4008	18	40	10	8	3832-4008	18	40	12	8
3810-4009	18	40	0	9	3816-4009	18	40	6	9	3830-4009	18	40	10	9	3832-4009	18	40	12	9
3810-4010	18	40	0	10	3816-4010	18	40	6	10	3830-4010	18	40	10	10	3832-4010	18	40	12	10
3810-4011	18	40	0	11	3816-4011	18	40	6	11	3830-4011	18	40	10	11	3832-4011	18	40	12	11
3810-4012	18	40	0	12	3816-4012	18	40	6	12	3830-4012	18	40	10	12	3832-4012	18	40	12	12
3810-4013	18	40	0	13	3816-4013	18	40	6	13	3830-4013	18	40	10	13	3832-4013	18	40	12	13
3810-4014	18	40	0	14	3816-4014	18	40	6	14	3830-4014	18	40	10	14	3832-4014	18	40	12	14
3810-4015	18	40	0	15	3816-4015	18	40	6	15	3830-4015	18	40	10	15	3832-4015	18	40	12	15
3810-4016	18	40	0	16	3816-4016	18	40	6	16	3830-4016	18	40	10	16	3832-4016	18	40	12	16
Part No.	Dimension (mm)				Part No.	Dimension (mm)				Part No.	Dimension (mm)				Part No.	Dimension (mm)			
	W(M-L)	L(A-P)	A(°)	H(mm)		W(M-L)	L(A-P)	A(°)	H(mm)		W(M-L)	L(A-P)	A(°)	H(mm)		W(M-L)	L(A-P)	A(°)	H(mm)
3810-4508	18	45	0	8	3816-4508	18	45	6	8	3830-4508	18	45	10	8	3832-4508	18	45	12	8
3810-4509	18	45	0	9	3816-4509	18	45	6	9	3830-4509	18	45	10	9	3832-4509	18	45	12	9
3810-4510	18	45	0	10	3816-4510	18	45	6	10	3830-4510	18	45	10	10	3832-4510	18	45	12	10
3810-4511	18	45	0	11	3816-4511	18	45	6	11	3830-4511	18	45	10	11	3832-4511	18	45	12	11
3810-4512	18	45	0	12	3816-4512	18	45	6	12	3830-4512	18	45	10	12	3832-4512	18	45	12	12
3810-4513	18	45	0	13	3816-4513	18	45	6	13	3830-4513	18	45	10	13	3832-4513	18	45	12	13
3810-4514	18	45	0	14	3816-4514	18	45	6	14	3830-4514	18	45	10	14	3832-4514	18	45	12	14
3810-4515	18	45	0	15	3816-4515	18	45	6	15	3830-4515	18	45	10	15	3832-4515	18	45	12	15
3810-4516	18	45	0	16	3816-4516	18	45	6	16	3830-4516	18	45	10	16	3832-4516	18	45	12	16
Part No.	Dimension (mm)				Part No.	Dimension (mm)				Part No.	Dimension (mm)				Part No.	Dimension (mm)			
	W(M-L)	L(A-P)	A(°)	H(mm)		W(M-L)	L(A-P)	A(°)	H(mm)		W(M-L)	L(A-P)	A(°)	H(mm)		W(M-L)	L(A-P)	A(°)	H(mm)
3810-5008	18	50	0	8	3816-5008	18	50	6	8	3830-5008	18	50	10	8	3832-5008	18	50	12	8
3810-5009	18	50	0	9	3816-5009	18	50	6	9	3830-5009	18	50	10	9	3832-5009	18	50	12	9
3810-5010	18	50	0	10	3816-5010	18	50	6	10	3830-5010	18	50	10	10	3832-5010	18	50	12	10
3810-5011	18	50	0	11	3816-5011	18	50	6	11	3830-5011	18	50	10	11	3832-5011	18	50	12	11
3810-5012	18	50	0	12	3816-5012	18	50	6	12	3830-5012	18	50	10	12	3832-5012	18	50	12	12
3810-5013	18	50	0	13	3816-5013	18	50	6	13	3830-5013	18	50	10	13	3832-5013	18	50	12	13
3810-5014	18	50	0	14	3816-5014	18	50	6	14	3830-5014	18	50	10	14	3832-5014	18	50	12	14
3810-5015	18	50	0	15	3816-5015	18	50	6	15	3830-5015	18	50	10	15	3832-5015	18	50	12	15
3810-5016	18	50	0	16	3816-5016	18	50	6	16	3830-5016	18	50	10	16	3832-5016	18	50	12	16
Part No.	Dimension (mm)				Part No.	Dimension (mm)				Part No.	Dimension (mm)				Part No.	Dimension (mm)			
	W(M-L)	L(A-P)	A(°)	H(mm)		W(M-L)	L(A-P)	A(°)	H(mm)		W(M-L)	L(A-P)	A(°)	H(mm)		W(M-L)	L(A-P)	A(°)	H(mm)
3810-5508	18	55	0	8	3816-5508	18	55	6	8	3830-5508	18	55	10	8	3832-5508	18	55	12	8
3810-5509	18	55	0	9	3816-5509	18	55	6	9	3830-5509	18	55	10	9	3832-5509	18	55	12	9
3810-5510	18	55	0	10	3816-5510	18	55	6	10	3830-5510	18	55	10	10	3832-5510	18	55	12	10
3810-5511	18	55	0	11	3816-5511	18	55	6	11	3830-5511	18	55	10	11	3832-5511	18	55	12	11
3810-5512	18	55	0	12	3816-5512	18	55	6	12	3830-5512	18	55	10	12	3832-5512	18	55	12	12
3810-5513	18	55	0	13	3816-5513	18	55	6	13	3830-5513	18	55	10	13	3832-5513	18	55	12	13
3810-5514	18	55	0	14	3816-5514	18	55	6	14	3830-5514	18	55	10	14	3832-5514	18	55	12	14
3810-5515	18	55	0	15	3816-5515	18	55	6	15	3830-5515	18	55	10	15	3832-5515	18	55	12	15
3810-5516	18	55	0	16	3816-5516	18	55	6	16	3830-5516	18	55	10	16	3832-5516	18	55	12	16
Part No.	Dimension (mm)				Part No.	Dimension (mm)				Part No.	Dimension (mm)				Part No.	Dimension (mm)			
	W(M-L)	L(A-P)	A(°)	H(mm)		W(M-L)	L(A-P)	A(°)	H(mm)		W(M-L)	L(A-P)	A(°)	H(mm)		W(M-L)	L(A-P)	A(°)	H(mm)
3810-6008	18	60	0	8	3816-6008	18	60	6	8	3830-6008	18	60	10	8	3832-6008	18	60	12	8
3810-6009	18	60	0	9	3816-6009	18	60	6	9	3830-6009	18	60	10	9	3832-6009	18	60	12	9
3810-6010	18	60	0	10	3816-6010	18	60	6	10	3830-6010	18	60	10	10	3832-6010	18	60	12	10
3810-6011	18	60	0	11	3816-6011	18	60	6	11	3830-6011	18	60	10	11	3832-6011	18	60	12	11
3810-6012	18	60	0	12	3816-6012	18	60	6	12	3830-6012	18	60	10	12	3832-6012	18	60	12	12
3810-6013	18	60	0	13	3816-6013	18	60	6	13	3830-6013	18	60	10	13	3832-6013	18	60	12	13
3810-6014	18	60	0	14	3816-6014	18	60	6	14	3830-6014	18	60	10	14	3832-6014	18	60	12	14
3810-6015	18	60	0	15	3816-6015	18	60	6	15	3830-6015	18	60	10	15	3832-6015	18	60	12	15
3810-6016	18	60	0	16	3816-6016	18	60	6	16	3830-6016	18	60	10	16	3832-6016	18	60	12	16

Instrument Specification

No.	Part No.	Description
1	LC04-0108	Trial 8 mm
2	LC04-0109	Trial 9 mm
3	LC04-0110	Trial 10 mm
4	LC04-0111	Trial 11 mm
5	LC04-0112	Trial 12 mm
6	LC04-0113	Trial 13 mm
7	LC04-0114	Trial 14 mm
8	LC04-0115	Trial 15 mm
9	LC04-0116	Trial 16 mm
10	LC04-0117	Trial 17 mm
11	LC04-0118	Trial 18 mm
12	LC04-0119	Trial 19 mm
13	LC04-0501	Distractor
14	LC04-1201	Curette Straight
15	LC04-1202	Curette Left Angled
16	LC04-1203	Curette Right Angled
17	LC04-1204	Curette Down Angled
18	LC04-1205	Curette Up Angled
19	LC04-0601	Chisel Straight
20	LC04-0602	Chisel Curved
21	LC04-0708	Broach 8x18 mm
22	LC04-0709	Broach 9x18 mm
23	LC04-0710	Broach 10x18 mm
24	LC04-0711	Broach 11x18 mm
25	LC04-0712	Broach 12x18 mm
26	LC04-0713	Broach 13x18 mm
27	LC04-0714	Broach 14x18 mm
28	LC04-0715	Broach 15x18 mm
29	LC04-0716	Broach 16x18 mm
30	LC04-0717	Broach 17x18 mm
31	LC04-0718	Broach 18x18 mm
32	LC04-0719	Broach 19x18 mm
33	LC04-0801	Insertor
34	LC04-0901	Implant Tamp
35	06TR0306	Modular T-Handle

Also available.



LnK TLIF Cage System



LnK PLIF Cage System

LnK MIS System



LnK Screw Fixation System

LnK Lumbar Intervertebral body Fusion Cage System

The LnK Lumbar Intervertebral body Fusion Cage System's implants are interbody fusion devices intended for use as an aid in spinal fixation. These hollow, rectangular implants are offered in a variety of widths, lengths, heights and lordotic angles designed to adapt to a variety of patient anatomies. They have serrations on the superior and inferior surfaces designed for fixation, ergonomically shaped anterior edges, and flat posterior edges. Radiopaque markers have been embedded within the implants, which are designed to allow for visualization in radiographic images

SURGICAL APPROACH

- PLIF(Posterior Lumbar Intervertebral body Fusion) PEEK Cages are to be implanted via posterior approach.
- TLIF(Transforaminal Lumbar Intervertebral body Fusion) PEEK Cages are to be implanted via transforaminal approach.
- ALIF(Anterior Lumbar Intervertebral body Fusion) PEEK Cages are to be implanted via anterior approach.

INDICATIONS

LnK Lumbar Intervertebral body Fusion Cage System is indicated for intervertebral body fusion procedures in skeletally mature patients with degenerative disc disease (DDD) at one or two contiguous levels from L2-S1. DDD is defined as discogenic back pain with degeneration of the disc confirmed by patient history and radiographic studies. These DDD patients may also have up to Grade 1 spondylolisthesis or retrolisthesis at the involved level(s). This device is to be used with autogenous bone graft. LnK Lumbar Intervertebral body Fusion Cage System is to be used with supplemental fixation. Patients should have at least six (6) months of non-operative treatment prior to treatment with an intervertebral cage.

GENERAL CONDITIONS OF USE

The implantation of intervertebral body fusion devices must be performed only by experienced spinal surgeons having undergone the necessary specific training in the use of such systems because this is a technically demanding procedure presenting a risk of serious injury to the patient.

CAUTION

Based on the fatigue testing results, the physician/surgeon must consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc. which may impact on the performance of the intervertebral body fusion device.

The implantation of the intervertebral body fusion device must be performed only by experienced spinal surgeons with specific training in the use of this device because this is a technically demanding procedure presenting a risk of serious injury to the patient.

Potential risks identified with the use of this intervertebral body fusion device, which may require additional surgery, include: device component fracture, loss of fixation, pseudoarthrosis (i.e. non-union), fracture of the vertebrae, neurological injury, and vascular or visceral injury.

Specialized instruments are provided by L&K Biomed and must be used to assure accurate implantation of the intervertebral body fusion device. While rare, intraoperative fracture or breakage of instruments can occur, instruments, which have experienced extensive use or extensive force, are more susceptible to fracture depending on the operative precaution, number of procedures, and disposal attention. Instruments must be examined for wear or damage prior to surgery. Instruments for implantation of the LnK Lumbar Intervertebral body Fusion Cage System is provided non-sterile and must be sterilized prior to use.

The LnK Lumbar Intervertebral body Fusion Cage System has not been evaluated for safety and compatibility in the MR environment. The LnK Lumbar Intervertebral body Fusion Cage System has not been tested for heating or migration in the MR environment.

INFECTION

Transient bacteremia can occur in daily life. Dental manipulation, endoscopic examination and other minor surgical procedures have been associated with transient bacteremia. To help prevent infection at the implant site, it may be advisable to use antibiotic prophylaxis before and after such procedures.

INSTRUMENTS

Specialized instruments are provided by L&K Biomed and must be used to assure accurate implantation of the device. While rare, intraoperative fracture or breakage of instruments can occur. Instruments which have experienced extensive use or extensive force are more susceptible to fracture depending on the operative precaution, number of procedures, disposal attention. Instruments must be examined for wear or damage prior to surgery.

REUSE

An implant should never be reused. While it may appear undamaged, a used implant may have acquired blemishes or latent compromise of its integrity which would reduce its service life.

Surgeons must verify that the instruments are in good condition and operating order prior to use during surgery.

HANDLING

Correct handling of the implant is extremely important. The operating surgeon must avoid notching or scratching the device.

ALLERGY AND HYPERSENSITIVITY TO FOREIGN BODIES

When hypersensitivity is suspected or proven, it is highly recommended that the tolerance of the skin to the materials that make up the implants be checked before they are implanted.

CONTRAINDICATIONS

Contraindications may be relative or absolute. The choice of a particular device must be carefully weighed against the patient's overall evaluation. Circumstances listed below may reduce the chances of a successful outcome:

- Prior fusion at the levels to be treated.
- An active infection at the operative site.
- Use except as indicated.
- Marked local inflammation.
- Any abnormality present which affects the normal process of bone remodeling including, but not limited to, severe osteoporosis involving the spine, bone absorption, osteopenia, primary or metastatic tumors involving the spine, active infection at the site or certain metabolic disorders affecting osteogenesis.
- Any mental or neuromuscular disorder which would create an unacceptable risk of fixation failure or complications in postoperative care.
- Open wounds.
- Pregnancy.
- Patients having inadequate tissue coverage of the operative site.
- Any neuromuscular deficit which places an unsafe load level on the device during the healing period.
- Obesity. An overweight or obese patient can produce loads on the spinal system which can lead to failure of the fixation of the device or to failure of the device itself. Obesity is defined according to the W.H.O. standards.
- A condition of senility, mental illness, or substance abuse. These conditions, among others, may cause the patient to ignore certain necessary limitations and precautions in the use of the implant, leading to failure or other complications.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests must be made prior to material selection or implantation.
- Other medical or surgical condition which would preclude the potential benefit of spinal implant surgery, such as the presence of tumors, congenital abnormalities, elevation of sedimentation rate unexplained by other diseases, elevation of white blood cell count (WBC), or marked left shift in the WBC differential count. These contra-indications can be relative or absolute and must be taken into account by the physician when making his decision. The above list is not exhaustive. Surgeons must discuss the relative contraindications with the patients.

INFORMATION FOR PATIENTS

The surgeon must discuss all physical and psychological limitations inherent to the use of the device with the patient. This includes the rehabilitation regimen, physical therapy, and wearing an appropriate orthosis as prescribed by the physician. Particular discussion should be directed to the issues of premature weight bearing, activity levels, and the necessity for periodic medical follow-up.

The surgeon must warn the patient of the surgical risks and made aware of possible adverse effects. The surgeon must warn the patient that the device cannot and does not replicate the flexibility, strength, reliability or durability of normal healthy bone, that the implant can break or become damaged as a result of strenuous activity or trauma, and that the device may need to be replaced in the future. If the patient is involved in an occupation or activity which applies inordinate stress upon the implant (e.g., substantial walking, running, lifting, or muscle strain) the surgeon must advise the patient that resultant forces can cause failure of the device. Patients who smoke have been shown to have an increased incidence of non-unions. Surgeons must advise patients of this fact and warn of the potential consequences. For diseased patients with degenerative disease, the progression of degenerative disease may be so advanced at the time of implantation that it may substantially decrease the expected useful life of the appliance. In such cases, orthopaedic devices may be considered only as a delaying technique or to provide temporary relief. Patients with previous spinal surgery at the level(s) to be treated may have different clinical outcomes compared to those without a previous surgery.

PREOPERATIVE PRECAUTIONS

- The surgical indication and the choice of implants must take into account certain important criteria such as:
- Patients involved in an occupation or activity that applies excessive loading upon the implant (e.g., substantial walking, running, lifting, or muscle strain) may be at increased risk for failure of the fusion and/or the device.
 - Surgeons must instruct patients in detail about the limitations of the implants, including, but not limited to, the impact of excessive loading through patient weight or activity, and be taught to govern their activities accordingly. The procedure will not restore function to the level expected with a normal, healthy spine, and the patient should not have unrealistic functional expectations.
 - A condition of senility, mental illness, chemical dependence or alcoholism. These conditions among others may cause the patients to ignore certain necessary limitations and precautions in the use of the implant, leading to failure and other complications.
 - Foreign body sensitivity. Where material sensitivity is suspected appropriate tests must be made prior to material implantation.
 - Surgeons must advise patients who smoke have been shown to have an increased incidence of non-unions. Such patients must be advised of this fact and warned of the potential consequences.
 - Care must be taken to protect the components from being marred, nicked, or notched as a result of contact with metal or abrasive objects.

THE CHOICE OF IMPLANTS

The choice of proper shape, size and design of the implant for each patient is crucial to the success of the surgery. The surgeon is responsible for this choice which depends on each patient. Patients who are overweight may be responsible for additional stresses and strains on the device which can speed up implant fatigue and/or lead to deformation or failure of the implants. The size and shape of the bone structures determine the size, shape and type of the implants. Once implanted, the implants are subjected to stresses and strains. These repeated stresses on the implants must be taken into consideration by the surgeon at the time of the choice of the implant, during implantation as well as in the post-operative follow-up period. Indeed, the stresses and strains on the implants may cause fatigue or fracture or deformation of the implants, before the bone graft has become completely consolidated. This may result in further side effects or necessitate the early removal of the osteosynthesis device.

INTRAOPERATIVE PRECAUTIONS

- The insertion of the implants must be carried out using instruments designed and provided for this purpose and in accordance with the specific implantation instructions for each implant. Those detailed instructions are provided in the surgical technique brochure supplied by L&K Biomed.
- Discard all damaged or mishandled implants.
- Never reuse an implant, even though it may appear undamaged.

POSTOPERATIVE PRECAUTIONS

Prior to adequate maturation of the fusion mass, implanted spinal instrumentation may need additional help to accommodate full load bearing. External support may be recommended by the physician from two to four months postoperatively or until x-rays or other procedures confirm adequate maturation of the fusion mass; external immobilization by bracing or casting be employed. Surgeons must instruct patients regarding appropriate and restricted activities during consolidation and maturation for the fusion mass in order to prevent placing excessive stress on the implants which may lead to fixation or implant failure and accompanying clinical problems. Surgeons must instruct patients to report any unusual changes of the operative site to his/her physician. The physician must closely monitor the patient if a change at the site has been detected. Patients with previous spinal surgery at the level(s) to be treated may have different clinical outcomes compared to those without a previous surgery.

ADVERSE EFFECTS

- Include but are not limited to:
 - Late bone fusion or no visible fusion mass and pseudarthrosis;
 - Peripheral neuropathies, nerve damage, heterotopic bone formation and neurovascular compromise, including paralysis;
 - While the expected life of spinal implant components is difficult to estimate, it is finite. These components are made of foreign materials which are placed within the body for the potential fusion of the spine and reduction of pain. However, due to the many biological, mechanical and physicochemical factors which affect these devices but cannot be evaluated in vivo, the components cannot be expected to indefinitely withstand the activity level and loads of normal healthy bone.
 - Superficial or deep-set infection and inflammatory phenomena;
 - Allergic reactions to the implanted materials although uncommon can occur;
 - Decrease in bone density due to stress shielding;
 - Dural leak requiring surgical repair.
 - Peripheral neuropathies, nerve damage, heterotopic bone formation and neurovascular compromise, including paralysis, loss of bowel or bladder function, or foot-drop may occur.
 - Cessation of growth of the fused portion of the spine.
 - Loss of proper spinal curvature, correction, height and/or reduction.
 - Delayed Union or Nonunion: Internal fixation appliances are load sharing devices which are used to obtain alignment until normal healing occurs. In the event that healing is delayed, does not occur, or failure to immobilize the delayed/nonunion results, the implant will be subject to excessive and repeated stresses which can eventually cause loosening, bending or fatigue fracture. The degree or success of union, loads produced by weight bearing, and activity levels will, among other conditions, dictate the longevity of the implant. If a nonunion develops or if the implants loosen, bend or break, the device(s) must be revised or removed immediately before serious injury occurs.
 - Neurological and spinal dura mater lesions from surgical trauma;
 - Early loosening may result from inadequate initial fixation, latent infection, premature loading of the device or trauma. Late loosening may result from trauma, infection, biological complications or mechanical problems, with the subsequent possibility of bone erosion, vertebral endplate injury or pain.
 - Serious complications may occur with any spinal surgery. These complications include, but are not limited to, genitourinary disorders; gastrointestinal disorders; vascular disorders, including thrombus; bronchopulmonary disorders, including emboli; bursitis, hemorrhage, myocardial infarction, infection, paralysis or death.
 - Inappropriate or improper surgical placement of this device may cause distraction or stress shielding of the graft or fusion mass. This may contribute to failure of an adequate fusion mass to form.
 - Intraoperative fissure, fracture, or perforation of the spine can occur due to implantation of bone graft or the intervertebral body above or below the level of surgery can occur due to trauma, the presence of defects, or poor bone stock. Adverse effects may necessitate reoperation. Adverse effects may necessitate reoperation or revision.
- The surgeon must warn the patient of these adverse effects as deemed necessary.

IMPLANT REMOVAL

- If fusion / bone graft growth occurs, the device will be deeply integrated into the bony tissues. As a result, the LnK Lumbar Intervertebral body Fusion Cage System is not intended to be removed unless the management of a complication or adverse event requires the removal. Any decision by a physician to remove the device should take into consideration such factors as:
- The risk to the patient of the additional surgical procedure as well as the difficulty of removal.
 - Migration of the implant, with subsequent pain and/or neurological, articular or soft tissue lesions.
 - Pain or abnormal sensations due to the presence of the implants.
 - Infection or inflammatory reactions.
 - Reduction in bone density due to the different distribution of mechanical and physiological stresses and strains.

CLEANING AND STERILIZATION

All implants and instruments used in surgery must be sterilized by the hospital prior to use. Remove all packaging materials prior to sterilization. Only sterile products should be placed in the operative field. Sterilization: recommended method to achieve a degree of sterility equal to at least 10⁻⁶. The gravity displacement sterilization parameters we suggested comply with AAMI ST79.L&K BIOMED recommends the following parameters:

METHOD	CYCLE	TEMPERATURE	EXPOSURE TIME
Steam	Gravity	270°F(132°C)	15Minutes (Dry time, 15-30 Minute)

Manufactured by:
L&K BIOMED Co.,Ltd.
1104-ho, 145, Gasandigital 1-ro, Seoul, 153-787 Korea
Tel. 82-2-2624-1471-4 / Fax. 82-2-2624-1477

SYMBOL TRANSLATION			
LOT NUMBER LOT	CATALOG NUMBER REF	DATE OF MANUFACTURE	SINGLE USE ONLY
NON-STERILE	MANUFACTURER	See package insert for labeling limitation	Federal Law (USA) restricts this device to sale, distribution, or use by or on the order of a physician... Rx Only



Distributed by Aegis Spine, Inc.
6400 S Fiddlers Green Circle Suite 1180, Greenwood Village, CO 80111, USA
Tel: +1-303-741-4123
Fax: +1-303-741-8131
Home page: <http://www.aegisspine.com>

Manufactured by LnK Biomed Co., Ltd.
#1104, ACE High-End Tower 3, 371-50, Gasan-dong,
Geumcheon-gu, Seoul, Korea
Tel : 82-2-2624-1471~4
Fax : 82-2-2624-1477

2013 LnK Biomed Korea, Inc.
All Rights Reserved.