




**Medtronic**

# CRESCENT<sup>®</sup> Spinal System Titanium Implants

Product Information



## Distinctive Design

- » Banana shape design
- » Bulleted nose
- » 6° Lordosis
- » Ball tip threaded inserter for inserting, tamping, and facilitating final positioning
- » Straight and 15° offset inserter attachment options
- » Teeth for fixation
- » 10mm and 12mm AP width for all implant lengths
- » Instrument set is applicable for MIS
- » Derotation tools to reposition an over-rotated compatible implant
- » Extraction tools to remove an implant



15° offset denoted by a "T" on the right side anterior portion of the implant

## Internal Volume

### Regular Width Implants

Part Number	Description	Total Internal Volume
6660725	7mm × 25mm	0.33cc
6660825	8mm × 25mm	0.37cc
6660925	9mm × 25mm	0.42cc
6661025	10mm × 25mm	0.46cc
6661125	11mm × 25mm	0.50cc
6661225	12mm × 25mm	0.55cc
6661325	13mm × 25mm	0.59cc
6661425	14mm × 25mm	0.63cc
6661525	15mm × 25mm	0.68cc
6661625	16mm × 25mm	0.72cc
<hr/>		
6660730	7mm × 30mm	0.50cc
6660830	8mm × 30mm	0.57cc
6660930	9mm × 30mm	0.68cc
6661030	10mm × 30mm	0.74cc
6661130	11mm × 30mm	0.89cc
6661230	12mm × 30mm	0.96cc
6661330	13mm × 30mm	1.11cc
6661430	14mm × 30mm	1.18cc
6661530	15mm × 30mm	1.33cc
6661630	16mm × 30mm	1.40cc
<hr/>		
6660736	7mm × 36mm	0.67cc
6660836	8mm × 36mm	0.76cc
6660936	9mm × 36mm	0.90cc
6661036	10mm × 36mm	0.99cc
6661136	11mm × 36mm	1.19cc
6661236	12mm × 36mm	1.28cc
6661336	13mm × 36mm	1.48cc
6661436	14mm × 36mm	1.57cc
6661536	15mm × 36mm	1.77cc
6661636	16mm × 36mm	1.87cc

### Wide Width Implants

Part Number	Description	Total Internal Volume
6670725	7mm × 25mm	0.50cc
6670825	8mm × 25mm	0.56cc
6670925	9mm × 25mm	0.62cc
6671025	10mm × 25mm	0.69cc
6671125	11mm × 25mm	0.75cc
6671225	12mm × 25mm	0.81cc
6671325	13mm × 25mm	0.87cc
6671425	14mm × 25mm	0.94cc
6671525	15mm × 25mm	1.00cc
6671625	16mm × 25mm	1.06cc
<hr/>		
6670730	7mm × 30mm	0.72cc
6670830	8mm × 30mm	0.81cc
6670930	9mm × 30mm	0.95cc
6671030	10mm × 30mm	1.05cc
6671130	11mm × 30mm	1.23cc
6671230	12mm × 30mm	1.33cc
6671330	13mm × 30mm	1.51cc
6671430	14mm × 30mm	1.61cc
6671530	15mm × 30mm	1.79cc
6671630	16mm × 30mm	1.89cc
<hr/>		
6670736	7mm × 36mm	1.04cc
6670836	8mm × 36mm	1.17cc
6670936	9mm × 36mm	1.38cc
6671036	10mm × 36mm	1.51cc
6671136	11mm × 36mm	1.80cc
6671236	12mm × 36mm	1.93cc
6671336	13mm × 36mm	2.21cc
6671436	14mm × 36mm	2.35cc
6671536	15mm × 36mm	2.63cc
6671636	16mm × 36mm	2.77cc

## Product Ordering Information

### SET TYPE 02347

#### CRESCENT® Spinal System - Titanium Implants

Part Number	Description	Set Quantity
6660725	7mm × 25mm	1
6660825	8mm × 25mm	2
6660925	9mm × 25mm	2
6661025	10mm × 25mm	2
6661125	11mm × 25mm	2
6661225	12mm × 25mm	2
6661325	13mm × 25mm	1
6661425	14mm × 25mm	1
6661525	15mm × 25mm	1
6661625	16mm × 25mm	1
6660730	7mm × 30mm	1
6660830	8mm × 30mm	1
6660930	9mm × 30mm	2
6661030	10mm × 30mm	2
6661130	11mm × 30mm	2
6661230	12mm × 30mm	2
6661330	13mm × 30mm	2
6661430	14mm × 30mm	1
6661530	15mm × 30mm	1
6661630	16mm × 30mm	1
6660736	7mm × 36mm	1
6660836	8mm × 36mm	1
6660936	9mm × 36mm	1
6661036	10mm × 36mm	1
6661136	11mm × 36mm	1
6661236	12mm × 36mm	1
6661336	13mm × 36mm	1
6661436	14mm × 36mm	1
6661536	15mm × 36mm	1
6661636	16mm × 36mm	1
	Suitcase	1

### SET TYPE 02348

#### CRESCENT® Spinal System – Titanium Implants, Wide

Part Number	Description	Set Quantity
6670725	7mm × 25mm	1
6670825	8mm × 25mm	2
6670925	9mm × 25mm	2
6671025	10mm × 25mm	2
6671125	11mm × 25mm	2
6671225	12mm × 25mm	2
6671325	13mm × 25mm	1
6671425	14mm × 25mm	1
6671525	15mm × 25mm	1
6671625	16mm × 25mm	1
6670730	7mm × 30mm	1
6670830	8mm × 30mm	1
6670930	9mm × 30mm	2
6671030	10mm × 30mm	2
6671130	11mm × 30mm	2
6671230	12mm × 30mm	2
6671330	13mm × 30mm	2
6671430	14mm × 30mm	1
6671530	15mm × 30mm	1
6671630	16mm × 30mm	1
6670736	7mm × 36mm	1
6670836	8mm × 36mm	1
6670936	9mm × 36mm	1
6671036	10mm × 36mm	1
6671136	11mm × 36mm	1
6671236	12mm × 36mm	1
6671336	13mm × 36mm	1
6671436	14mm × 36mm	1
6671536	15mm × 36mm	1
6671636	16mm × 36mm	1
	Suitcase	1

[www.medtronic.com](http://www.medtronic.com)

**Medtronic**

Spinal and Biologics Business  
Worldwide Headquarters

2600 Sofamor Danek Drive  
Memphis, TN 38132

1801 Pyramid Place  
Memphis, TN 38132

(901) 396-3133  
(800) 876-3133  
Customer Service: (800) 933-2635

For more information visit  
[www.myspinetools.com](http://www.myspinetools.com)

**Medtronic Europe Headquarters**

Medtronic International Trading SÀRL

Route de Molliau 31  
1131 Tolochenaz  
Switzerland

Tel: +41-21-802-7000  
Fax: +41-21-802-7275

Please see the package insert for the complete list of indications, warnings, precautions, and other important medical information.

Potential risks associated with the device include, but are not limited to:

- » Implant migration
- » Loss of spinal curvature, correction, height, and/or reduction
- » Bone fracture or stress shielding at, above, or below the level of surgery
- » Bone graft donor site complication
- » Loss of, or increase in, spinal mobility or function



**Medtronic**