



ODALYS

Posterior thoraco-lumbar
osteosynthesis system





List of **ODALYS** instruments




Ratchet Handle - **SP0040** - (x2)



Torque Control T-handle - **SP0030** - (x1*)



Rod Gripper - **PC0040** - (x1)



Axial Persuader - **GU0030** - (x1)



Counter-torque - **GU0020** - (x1)



Transverse connector torque control handle - **SP0060** - (x1*)



Screw Holder - **SD0100** - (x2)



Short Holding Screwdriver T25 - **SD0050** - (x1)



Long Holding Screwdriver T25 - **SD0020** - (x2)



Thin Straight Feeler - **PA0010** - (x1)



Thin curved feeler - **PA0020** - (x1)



Malleable Rod L60 - **GB0050** - (x1)



Malleable Rod L120 - **GB0060** - (x1)



Pedicular Marker - **PA0040** - (x4)



Dual Ball-ended Pedicular Marker - **PA0050** - (x4)



Rod Holder - **PC0030** - (x1)



De-rotation Wrench - **WR0010** - (x1)



Curved Thoracic Probe - **PA0030** - (x1)



Long Round Awl - **PN0050** - (x1)



Rocker - **DV0060** - (x1)



Screw Tap Ø4.2mm - **TD0040** - (x1)



Screw Tap Ø5mm - **TD0010** - (x1)



Screw Tap Ø6mm - **TD0020** - (x1)



Screw Tap Ø7mm - **TD0030** - (x1)



Compressor - **PC0050** - (x1)



French Bender - **DV0130** - (x1)



Distractor - **PC0060** - (x1)

* Instrument with limited shelf life end. Please consult the provided inventory or contact KISCO International to the date of guarantee.



List of **ODALYS** optional instruments



Short Screwdriver T25 - [SD0060](#)



Straight Thoracic Probe - [PA0070](#)



Short Square Awl - [PN0040](#)



Variable Square Awl - [PN0060](#)




Rod Pusher - [DV0080](#)



Screw Head Positionner - [DV0090](#)



Instrument/implant box lid - BT0310 - (x2)



Instrument box 1 - BT0380 - (x1)



Instrument box 2 - BT0390 - (x1)



Instrument box 3 - BT0400 - (x1)

Dual ball-ended Pedicular Marker (Stainless Steel) - PA0050 - {x4}

Curved Thoracic Probe (Stainless Steel + Silicone) - PA0030 - {x1}
[option : PA0070]


Pedicular Marker (Stainless Steel) - PA0040 - {x4}


Long Round Awl
(Stainless Steel + Silicone)
PN0050 - {x1}
[option : PN0040 or
PN0060]

Ratchet Handle (Stainless Steel + Silicone) - SP0040 - {x2}


Screw Holder (Stainless Steel + Radel) - SD0100 - {x2}

Instrument Box 1 - BT0380 - (x1)







Intrument Box 2 - BT0390 - (x1)



* Instrument with limited shelf life end. Please consult the provided inventory or contact KISCO International to the date of guarantee.







List of **ODALYS** sterile implants Ø6

2 x Polyaxial Screws Ø5

L30	L35	L40	L45
[x2]	[x3]	[x3]	[x2]
10PS0010	10PS0020	10PS0030	10PS0040



2 x Polyaxial Screws Ø6

L30	L35	L40	L45	L50
[x2]	[x4]	[x4]	[x4]	[x3]
10PS0050	10PS0060	10PS0070	10PS0080	10PS0090



2 x Polyaxial Screws Ø7

L30	L35	L40	L45	L50
[x1]	[x2]	[x2]	[x2]	[x1]
10PS0100	10PS0110	10PS0120	10PS0130	10PS0140



2 x Polyaxial Spendylolisthesis Screws Ø6

L30	L35	L40	L45	L50
[x1]	[x1]	[x2]	[x1]	[x1]
10SS0010	10SS0020	10SS0030	10SS0040	10SS0050



2 x Polyaxial Spendylolisthesis Screws Ø7

L30	L35	L40	L45	L50
[x1]	[x1]	[x1]	[x1]	[x1]
10SS0060	10SS0070	10SS0080	10SS0090	10SS0100

2 x Monoaxial Screws Ø4.2

L25	L30	L35	L40
[x2]	[x2]	[x2]	[x2]
10MS0010	10MS0020	10MS0030	10MS0040



2 x Monoaxial Screws Ø5

L30	L35	L40	L45
[x1]	[x1]	[x2]	[x1]
10MS0050	10MS0060	10MS0070	10MS0080



2 x Monoaxial Screws Ø6

L35	L40	L45
[x1]	[x3]	[x2]
10MS0160	10MS0170	10MS0180





2 x Curved Rods Ø6

L35	L40	L45	L50	L55	L65	L75	L85	L95	L105
[x2]	[x2]	[x1]	[x1]	[x1]	[x1]	[x1]	[x1]	[x1]	[x1]
10BR0010	10BR0020	10BR0030	10BR0040	10BR0050	10BR0060	10BR0070	10BR0080	10BR0090	10BR0100



Transverse system

L37-49	L49-73
[x2]	[x2]
10CLO010	10CLO020



2 x Locking Screws

[x3]
10NT0010





List of **ODALYS** non-sterile implants Ø5.5

Polyaxial Screw Ø4

L25	L30	L35	L40	L45	L50
x4	x4	x4	*	*	*
ODPS55425	ODPS55430	ODPS55435	ODPS55440	ODPS55445	ODPS55450



Polyaxial Screw Ø5

L25	L30	L35	L40	L45	L50	L55
*	x4	x6	x6	x4	x2	*
ODPS55525	ODPS55530	ODPS55535	ODPS55540	ODPS55545	ODPS55550	ODPS55555



Polyaxial Screw Ø6

L25	L30	L35	L40	L45	L50	L55
*	x2	x6	x6	x6	x6	x6
ODPS55625	ODPS55630	ODPS55635	ODPS55640	ODPS55645	ODPS55650	ODPS55655



Polyaxial Screw Ø7

L25	L30	L35	L40	L45	L50	L55
*	x2	x4	x2	x4	x4	x2
ODPS55725	ODPS55730	ODPS55735	ODPS55740	ODPS55745	ODPS55750	ODPS55755



Monoaxial Screw Ø4.2

L25	L30	L35	L40	L45	L50
x2	x4	x4	*	*	*
ODMS55425	ODMS55430	ODMS55435	ODMS55440	ODMS55445	ODMS55450



Monoaxial Screw Ø5

L25	L30	L35	L40	L45	L50	L55
*	x2	x4	x6	x4	x2	*
ODMS55525	ODMS55530	ODMS55535	ODMS55540	ODMS55545	ODMS55550	ODMS55555



Monoaxial Screw Ø6

L25	L30	L35	L40	L45	L50	L55
*	*	x2	x6	x6	x2	x2
ODMS55625	ODMS55630	ODMS55635	ODMS55640	ODMS55645	ODMS55650	ODMS55655



Monoaxial Screw Ø7

L25	L30	L35	L40	L45	L50	L55
*	*	*	*	x2	x4	x2
ODMS55725	ODMS55730	ODMS55735	ODMS55740	ODMS55745	ODMS55750	ODMS55755

* On request, conditions apply



List of **ODALYS** non-sterile implants Ø5.5

Polyaxial Spondylolisthesis Screw Ø4

L25	L30	L35	L40	L45	L50
x2	x4	x4	x2	*	*
ODSS55425	ODSS55430	ODSS55435	ODSS55440	ODSS55445	ODSS55450



Polyaxial Spondylolisthesis Screw Ø5

L25	L30	L35	L40	L45	L50	L55
*	x2	x4	x6	x6	x4	*
ODSS55525	ODSS55530	ODSS55535	ODSS55540	ODSS55545	ODSS55550	ODSS55555



Polyaxial Spondylolisthesis Screw Ø6

L25	L30	L35	L40	L45	L50	L55
*	x2	x4	x6	x6	x6	x4
ODSS55625	ODSS55630	ODSS55635	ODSS55640	ODSS55645	ODSS55650	ODSS55655



Polyaxial Spondylolisthesis Screw Ø7

L25	L30	L35	L40	L45	L50	L55
*	*	x2	x2	x4	x4	x2
ODSS55725	ODSS55730	ODSS55735	ODSS55740	ODSS55745	ODSS55750	ODSS55755

Monoaxial Spondylolisthesis Screw Ø4.2

L25	L30	L35	L40	L45	L50
*	*	*	*	*	*
ODSMS55425	ODSMS55430	ODSMS55435	ODSMS55440	ODSMS55445	ODSMS55450



Monoaxial Spondylolisthesis Screw Ø5

L25	L30	L35	L40	L45	L50
*	*	*	*	*	*
ODSMS55525	ODSMS55530	ODSMS55535	ODSMS55540	ODSMS55545	ODSMS55550

Monoaxial Spondylolisthesis Screw Ø6

L25	L30	L35	L40	L45	L50	L55
*	*	*	*	*	*	*
ODSMS55625	ODSMS55630	ODSMS55635	ODSMS55640	ODSMS55645	ODSMS55650	ODSMS55655



Monoaxial Spondylolisthesis Screw Ø7

L25	L30	L35	L40	L45	L50	L55
*	*	*	*	*	*	*
ODSMS55725	ODSMS55730	ODSMS55735	ODSMS55740	ODSMS55745	ODSMS55750	ODSMS55755

* On request, conditions apply



List of **ODALYS** non-sterile implants Ø5.5

Curved Rod Ø5.5

L35	L45	L55	L65	L75	L85	L95	L105	L115
x2	x2	x2	x2	x2	x2	x2	x2	x2
ODPBR55035	ODPBR55045	ODPBR55055	ODPBR55065	ODPBR55075	ODPBR55085	ODPBR55095	ODPBR55105	ODPBR55115



Rod Ø5.5

L135	L155	L185	L225	L255
x2	x2	x2	x2	x2
ODRR55135	ODRR55155	ODRR55185	ODRR55225	ODRR55255

Transverse system


L31-40	L40-57	L57-91
x2	x2	x2
ODTLD550	ODTLD551	ODTLD552



Locking Cap

x30
ODLS







List of **ODALYS** non-sterile implants Ø6

Polyaxial Screw Ø4

L25	L30	L35	L40	L45	L50
x4	x4	x4	*	*	*
ODPS6425	ODPS6430	ODPS6435	ODPS6440	ODPS6445	ODPS6450



Polyaxial Screw Ø5

L25	L30	L35	L40	L45	L50	L55
*	x4	x6	x6	x4	x2	*
ODPS6525	ODPS6530	ODPS6535	ODPS6540	ODPS6545	ODPS6550	ODPS6555

Polyaxial Screw Ø6

L25	L30	L35	L40	L45	L50	L55
*	x2	x6	x6	x6	x6	x6
ODPS6625	ODPS6630	ODPS6635	ODPS6640	ODPS6645	ODPS6650	ODPS6655



Polyaxial Screw Ø7

L25	L30	L35	L40	L45	L50	L55
*	x2	x4	x2	x4	x4	x2
ODPS6725	ODPS6730	ODPS6735	ODPS6740	ODPS6745	ODPS6750	ODPS6755

Monoaxial Screw Ø4.2

L25	L30	L35	L40	L45	L50
x2	x4	x4	*	*	*
ODMS6425	ODMS6430	ODMS6435	ODMS6440	ODMS6445	ODMS6450



Monoaxial Screw Ø5

L25	L30	L35	L40	L45	L50	L55
*	x2	x4	x6	x4	x2	*
ODMS6525	ODMS6530	ODMS6535	ODMS6540	ODMS6545	ODMS6550	ODMS6555

Monoaxial Screw Ø6

L25	L30	L35	L40	L45	L50	L55
*	*	x2	x6	x6	x2	x2
ODMS6625	ODMS6630	ODMS6635	ODMS6640	ODMS6645	ODMS6650	ODMS6655



Monoaxial Screw Ø7

L25	L30	L35	L40	L45	L50	L55
*	*	*	*	x2	x4	x2
ODMS6725	ODMS6730	ODMS6735	ODMS6740	ODMS6745	ODMS6750	ODMS6755

* On request, conditions apply



List of **ODALYS** non-sterile implants Ø6

Polyaxial Spondylolisthesis Screw Ø4

L25	L30	L35	L40	L45	L50
x2	x4	x4	x2	*	*
ODSS6425	ODSS6430	ODSS6435	ODSS6440	ODSS6445	ODSS6450



Polyaxial Spondylolisthesis Screw Ø5

L25	L30	L35	L40	L45	L50	L55
*	x2	x4	x6	x6	x4	*
ODSS6525	ODSS6530	ODSS6535	ODSS6540	ODSS6545	ODSS6550	ODSS6555

Polyaxial Spondylolisthesis Screw Ø6

L25	L30	L35	L40	L45	L50	L55
*	x2	x4	x6	x6	x6	x4
ODSS6625	ODSS6630	ODSS6635	ODSS6640	ODSS6645	ODSS6650	ODSS6655



Polyaxial Spondylolisthesis Screw Ø7

L25	L30	L35	L40	L45	L50	L55
*	*	x2	x2	x4	x4	x2
ODSS6725	ODSS6730	ODSS6735	ODSS6740	ODSS6745	ODSS6750	ODSS6755

Monoaxial Spondylolisthesis Screw Ø4.2

L25	L30	L35	L40	L45	L50
*	*	*	*	*	*
ODSMS6425	ODSMS6430	ODSMS6435	ODSMS6440	ODSMS6445	ODSMS6450



Monoaxial Spondylolisthesis Screw Ø5

L25	L30	L35	L40	L45	L50
*	*	*	*	*	*
ODSMS6525	ODSMS6530	ODSMS6535	ODSMS6540	ODSMS6545	ODSMS6550

Monoaxial Spondylolisthesis Screw Ø6

L25	L30	L35	L40	L45	L50	L55
*	*	*	*	*	*	*
ODSMS6625	ODSMS6630	ODSMS6635	ODSMS6640	ODSMS6645	ODSMS6650	ODSMS6655



Monoaxial Spondylolisthesis Screw Ø7

L25	L30	L35	L40	L45	L50	L55
*	*	*	*	*	*	*
ODSMS6725	ODSMS6730	ODSMS6735	ODSMS6740	ODSMS6745	ODSMS6750	ODSMS6755

* On request, conditions apply



List of **ODALYS** non-sterile implants Ø6

Curved Rod Ø6

L35	L45	L55	L65	L75	L85	L95	L105	L115
x2	x2	x2	x2	x2	x2	x2	x2	x2
ODPBR6035	ODPBR6045	ODPBR6055	ODPBR6065	ODPBR6075	ODPBR6085	ODPBR6095	ODPBR6105	ODPBR6115



Rod Ø6

L135	L155	L185	L225	L255
x2	x2	x2	x2	x2
ODRR6135	ODRR6155	ODRR6185	ODRR6225	ODRR6255

Transverse System


L31-40	L40-57	L57-91
x2	x2	x2
ODTLD60	ODTLD61	ODTLD62



Locking Cap

x30
ODLS








List of **ODALYS** non-sterile implant box



Implant Box - **BT0300** - {x1}

Instrument/Implant Box Lid - **BT0310** - {x1}



Screw Rack - **BT0320** - {x1}

Reduction* Screw Rack - **BT0330** - {x1}



* Spondylolisthesis



List of **ODALYS** non-sterile implant box

Locking Cap Box - **BT0340** - (x1)



Locking cap box tray - **BT0345** - (x1)

Locking Cap Box Lid - **BT0350** - (x1)



Transverse Connector Box - **BT0360** - (x1)

Transverse Connector Box Lid - **BT0370** - (x1)





List of ODALYS DEFORMITY non-sterile implants

Pedicle hook Ø5.5

Small	Medium	Large
x6	x6	x6
10HK0110	10HK0111	10HK0112

Pedicle hook Ø6

Small	Medium	Large
x6	x6	x6
10HK0115	10HK0116	10HK0117



Laminar hook Ø5.5

Small	Large
x6	x6
10HK0120	10HK0121

Laminar hook Ø6

Small	Large
x6	x6
10HK0125	10HK0126

Extended hook Ø5.5

Small	Large
x4	x4
10HK0130	10HK0131

Extended hook Ø6

Small	Large
x4	x4
10HK0135	10HK0136



Angled hook Ø5.5

Small Left	Large Left	Small Right	Large Right
x4	x4	x4	x4
10HK0140	10HK0141	10HK0142	10HK0143

Angled hook Ø6

Small Left	Large Left	Small Right	Large Right
x4	x4	x4	x4
10HK0145	10HK0146	10HK0147	10HK0148

Offset hook Ø5.5

Small Left	Large Left	Small Right	Large Right
x2	x2	x2	x2
10HK0150	10HK0151	10HK0152	10HK0153

Offset hook Ø6

Small Left	Large Left	Small Right	Large Right
x2	x2	x2	x2
10HK0155	10HK0156	10HK0157	10HK0158





List of **ODALYS DEFORMITY** non-sterile implants

Rod Ø5.5

L285	L345	L385	L425	L505
x2	x2	x2	x2	x2
ODRR55285	ODRR55345	ODRR55385	ODRR55425	ODRR55505



Rod Ø6

L285	L345	L385	L425	L505
x2	x2	x2	x2	x2
ODRR6285	ODRR6345	ODRR6385	ODRR6425	ODRR6505

Locking Cap

x30

ODLS





List of ODALYS DEFORMITY instruments



Rod gripper - PC0040 - (x1)



Coronal bender Ø5.5
left : FC0025 - (x1)
right : FC0020 - (x1)

OR

Coronal bender Ø6
left : FC0045 - (x1)
right : FC0040 - (x1)



In situ bender Ø5.5
left - FC0015 - (x1)
right - FC0010 - (x1)

OR

In situ bender Ø6
left - FC0035 - (x1)
right - FC0030 - (x1)



Pedicular starter - RG0010 - (x1)



Transverse process starter - RG0020 - (x1)



Laminar starter - RG0030 - (x1)



Hook impactor - AJ0030 - (x1)



Hook pusher - DV0120 - (x1)



Malleable rod L500 - GB0080 - (x1)



Straight hook holder - HD0070 - (x1)



Curved hook holder - HD0071 - (x1)




Pedicle hook Ø5.5 (Titanium)

Small	Medium	Large
x6	x6	x6
10HK0110	10HK0111	10HK0112

OR

Pedicle hook Ø6 (Titanium)

Small	Medium	Large
x6	x6	x6
10HK0115	10HK0116	10HK0117



Laminar hook Ø5.5 (Titanium)

Small	Large
x6	x6
10HK0120	10HK0121

OR

Laminar hook Ø6 (Titanium)

Small	Large
x6	x6
10HK0125	10HK0126

Extended hook Ø5.5 (Titanium)


Small	Large
x4	x4
10HK0130	10HK0131

OR

Extended hook Ø6 (Titanium)

Small	Large
x4	x4
10HK0135	10HK0136

Rod rack - BT0610 - (x2)
 Hook rack - BT0620 - (x1) and BT0630 - (x1)



Angled hook Ø5.5 (Titanium)

Small Left	Large Left	Small Right	Large Right
x4	x4	x4	x4
10HK0140	10HK0141	10HK0142	10HK0143

OR

Angled hook Ø6 (Titanium)

Small Left	Large Left	Small Right	Large Right
x4	x4	x4	x4
10HK0145	10HK0146	10HK0147	10HK0148

Rod Ø5.5 (Titanium)

L285	L345	L385	L425	L505
x2	x2	x2	x2	x2
ODRR55285	ODRR55345	ODRR55385	ODRR55425	ODRR55505

OR

Rod Ø6 (Titanium)

L285	L345	L385	L425	L505
x2	x2	x2	x2	x2
ODRR6285	ODRR6345	ODRR6385	ODRR6425	ODRR6505


Offset hook Ø5.5 (Titanium)

Small Left	Large Left	Small Right	Large Right
x2	x2	x2	x2
10HK0150	10HK0151	10HK0152	10HK0153

OR


Offset hook Ø6 (Titanium)

Small Left	Large Left	Small Right	Large Right
x2	x2	x2	x2
10HK0155	10HK0156	10HK0157	10HK0158



Deformity instrument/implant box - BT0600 - (x1)

Transverse process starter (Stainless Steel + Silicone) - RG0020 - (x1)



In situ bender Ø5.5
(Stainless Steel)
left - FC0015 - (x1)
right - FC0010 - (x1)

Hook pusher
(Stainless Steel + Silicone)
DVD120 - (x1)

Laminar starter
(Stainless Steel + Silicone)
RG0030 - (x1)

Hook impactor - (Stainless Steel + Silicone) - AJ0030 - (x1)

Coronal bender Ø5.5
(Stainless Steel)
left : FC0025 - (x1)
right : FC0020 - (x1)

OR

Coronal bender Ø6
(Stainless Steel)
left : FC0045 - (x1)
right : FC0040 - (x1)

In situ bender Ø6 (Stainless Steel)
left - FC0035 - (x1)
right - FC0030 - (x1)

