

Portico[™]

TLIF CAGE SYSTEM




System Features

- Angled Vents Allow for a Controlled Flow of Graft Material From the Device Into the Anterior Disc Space
- Bulleted Tip Provides Ease of Insertion and Precise Placement
- Aggressive Tooth Design Prevents Backing Out
- Tantalum Rod Markers Provide Excellent Visualization Under Fluoroscopy for Precise Placement
- Large Load Bearing Surface Area Offers Greater Contact with the Apophyseal Ring, Preventing Subsidence



Crafted for Perfection Portico™

TLIF Cage



Height



Footprint



Portico Convex TLIF

		Length (mm)	
		28	Graft Volume (cc)
Height (mm)	8	✓	0.75
	9	✓	0.85
	10	✓	0.95
	11	✓	1.0
	12	✓	1.1
	14	✓	1.3

■ Standard Sizes

With **InFill** Technology
Interbody Fusion System